

CEL MXD

Independence Solution

Tunstall

What is it?

The Tunstall CEL MXD is a Medical Alert System which utilizes the AT&T cellular network to communicate with our Emergency Response Center. A landline or traditional phone service is not required, nor does it require a cellular plan. The CEL MXD system includes a two-way voice base console and one personal activator which can be worn as a lanyard, or wristband.

How does it work?

Subscribers call for help by pressing the button on their personal activator or CEL MXD console. The call is sent to our 24-hour Emergency Response Center where trained operators can access call details (subscriber's name, address, and relevant information) that are displayed on their screen.

Our Emergency Response Center operator communicates with the subscriber through the CEL MXD console's loudspeaker and microphone to determine the nature of the call and the appropriate course of action.



Who is it for?

A wide variety of subscribers can benefit from the CEL MXD's state-of-the-art technology, including people who:

- Have been discharged from the hospital and require additional support and assistance at home
- Are of any age living alone
- Are older or have limited mobility
- Are subject to domestic violence or harassment
- Have household security issues or are living in fear of crime
- Are caring for someone in their home
- Have long-term conditions such as dementia, heart disease, or COPD, etc.
- Have mental, physical, or sensory impairments

The CEL MXD offers family, caregivers, and friends the reassurance that their loved ones can use their personal activator to easily call for help whenever they need it in and around their home.

CEL MXD

Independence Solution

Key Features

- Professional Emergency Response Centers
- Audible announcements available in English or Spanish
- 60 hour back-up console battery
- Integrated loud speaker & microphone
- Operates on the AT&T cellular network*
- Supports multiple personal activators
- FDA and ETL listed
- Compatible with the MyActive Fall Detector**

Benefits

- Available 24/7, 365
- Subscriber is informed of the call progress and status
- Provides reassurance in the event of a power outage
- Enables two-way voice communication between subscriber and Emergency Call Center Operator
- No need for land-line, no cellular plan needed. Operates on a leading reliable network.
- Allows multiple subscribers to use a single console unit
- Tunstall CEL MXD is FDA registered, complies with all medical device regulations, and is ETL listed to conform with UL 1635 and UL 1637 product safety standards
- Automatically senses when a fall occurs

Personal Activator

- Multiple wearing options
- 600ft activator range
- Water resistant (IP67)
- Necklace or bracelet
- Range suitable for typical dwellings/homes
- Can be worn in shower or bath

Specifications

- | | |
|-----------------------------|---|
| • Part(s) for this solution | • CEL MXD (MT100) Standard Activator: (MT004)
MyActive Fall Detector: (MT002) |
| • Dimensions | • 2.5" x 5.25" x 7.8" (HxWxL) |
| • Weight | • 1.8 lbs |
| • Battery | • 60 hour back-up base battery |
| • Installation Requirements | • AC power outlet
AT&T 3G cellular coverage* |
| • Maintenance Requirements | • Periodic testing |

For additional information please call

1.800.286.2622

or email US.Info@tunstall.com

36-36 33rd Street, Suite 103, Long Island City, NY 11106

americas.tunstall.com

*The Tunstall CEL MXD requires a satisfactory AT&T Cellular network connection to operate as intended. Coverage may not be available in all areas and may change after system activation due to a number of factors.

**Fall detection feature does not detect 100% of falls. If able, users should always press their help button when they need help.

MA-16-0225S Rev 2.0

Tunstall