

# belle®

Freedom Solution

Tunstall

SOLUTION SHEET

## What is it?

The Tunstall belle® is a small, lightweight, mobile Medical Alert System that enables seniors and individuals with special needs to live independently and call for help anytime, anywhere in the U.S. where there is cellular coverage. The product features simple one button access, two-way voice communication, Location Assistance Technology,\* and a battery that lasts up to 30 days on a single charge.\*\*

## How does it work?

Subscribers call for help, anywhere or anytime, by pressing the belle® call button. The call is sent to our 24-hour Emergency Response Center (ERC) where trained operators can access call details (subscriber's name, address, and relevant information) that are displayed on their screen.

If a subscriber is unable to confirm their location during an emergency, the device's location assistance technology can be used by dispatchers to help direct responders to the subscriber.



## Who is it for?

A wide variety of subscribers can benefit from the features of the belle®, including people who:

- Feel unsafe outside the home
- Have household security issues or are living in fear of crime
- Want to lead an active lifestyle
- Are subject to domestic violence or harassment
- Are of any age living alone
- Travel alone

The belle® gives subscribers the confidence they need when venturing outside the home; providing peace of mind and independence.

## Key Features Benefits

- |   |  |
|---|--|
| <ul style="list-style-type: none"> <li>• Simple, one-button access</li> <li>• Convenient user friendly charging cradle</li> <li>• AT&amp;T® or Verizon® cellular network***</li> <li>• Rechargeable 30 day standby battery**</li> <li>• Integrated loud speaker &amp; microphone</li> <li>• Location Assistance Technology</li> <li>• 911 back-up calling</li> <li>• Water resistant (IP67)</li> <li>• Compact, lightweight design</li> </ul> | <ul style="list-style-type: none"> <li>• Quickly connects subscribers to our trained Emergency Response Center Operators</li> <li>• Simple design lets subscribers know the device is properly docked and can be located next to the bed for charging and easy access to help at night</li> <li>• AT&amp;T and Verizon are leading cellular network providers with enhanced reliability. Not affected by the 2G network sunset.</li> <li>• Extended standby battery life affords subscribers the luxury of charging the device when convenient and reduces the likelihood of a dead battery</li> <li>• Enables two-way voice communication between subscriber and Emergency Call Center Operator</li> <li>• Cellular and wifi network based mobile technologies which work indoors or outdoors and can assist our Operators and 911 dispatchers in the location of subscribers</li> <li>• If the device cannot connect via the AT&amp;T® or Verizon® network it will attempt to call 911 directly on an alternate compatible carrier network***</li> <li>• Can be worn in shower or bath</li> <li>• Easy to handle and can be discreetly worn or carried in a pocket or purse</li> </ul> |
|---|--|

## Specifications

- |   |  |
|---|--|
| <ul style="list-style-type: none"> <li>• Part(s) for this solution</li> <li>• Dimensions</li> <li>• Weight</li> <li>• Battery</li> <li>• Installation Requirements</li> <li>• Maintenance Requirements</li> </ul> | <ul style="list-style-type: none"> <li>• belle® (196-001) AT&amp;T (196-101) Verizon</li> <li>• 2.7" x 1.7" x 0.7" (WxLxD)</li> <li>• 2.0 oz</li> <li>• 30 days standby battery (rechargeable)</li> <li>• AC power outlet<br/>AT&amp;T® or Verizon® cellular coverage***</li> <li>• Periodic charging of the device<br/>Gentle cleaning</li> </ul> |
|---|--|

For additional information please call  
**1.800.286.2622**  
 or email [USInfo@tunstall.com](mailto:USInfo@tunstall.com)  
 36-36 33rd Street, Suite 103, Long Island City, NY 11106  
[americas.tunstall.com](http://americas.tunstall.com)

\*Location accuracy may vary depending on network conditions and availability  
 \*\*30 day battery life based on standard operation and may vary with activity  
 \*\*\* AT&T and Verizon are trademarks of AT&T and Verizon  
 belle® is a Registered trademark of FreeUS

